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
- 界面会弹出以下信息(开机后, 滴定头不在起始位置, 请确定执行)
- The buret has not been Home Reset since the instrument was powered, Please click ok to Home Reset the buret.

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
Untitled - ITCRun

File View Experiment Buret Help Main Menu 主菜单

Buttons 快捷键

25.0000 C -159.4811 μ W 58772 s

Setup Monitor Data System Diagnostics Page Tabs 页面选择

Nano ITC

Run Mode Indicators 运行模式指示

ONLINE EQUILIBRATE PROGRAM

Stirring rate (RPM): 250 ☐ Stirrer on 搅拌转速

Syringe Size (μ L): 50

Experiment Type

☒ Incremental Titration

☐ Continuous Titration

☐ Electrical Pulses

Temperature Control

Set Point: 25

Apply

Data interval (s): 1

Data protection

☒ Save Experiment

Interval (minutes): 1

样品数据 Sample Data

Experiment Details

Syringe Concentration (mM): 1.3

Cell Concentration (mM): 0.07

Comments: RNase
2-CMP
250 rpm
25 degrees C

Instrument: Nano ITC 2G

Firmware Version: 1.18

Experiment Event Timing and Equilibration Settings 实验步骤及平衡设置

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Save

Load

☒ Start Delay (s): 0 (Equilibration time without collecting data)

☐ Auto Equilibrate

Expected Heats:

☒ Small ☐ Med ☐ Large

☐ Timeout (s): 0

Initial Baseline (s): 250 (Baseline collected before first injection)

Final Baseline (s): 0 (Baseline collected after last injection)

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The screenshot displays the 'Setup' tab of the Nano-ITC control software. It features several sections for configuring the experiment:

- Instrument Control:** Contains fields for 'Stirring Rate (RPM)' (set to 250), 'Syringe Size (μL)' (set to 250), 'Temperature Set Point (°C)' (set to 25), and 'Data Interval' (set to 5). Each field has an associated 'Update' button with a circular arrow icon.
- Auto Save Experiment:** Includes a checked checkbox for 'Save data every' followed by a field set to '10' minutes.
- Experiment Details:** Contains input fields for 'Syringe Concentration (mM)' and 'Cell Concentration (mM)', and a 'Comments' text area.

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
The screenshot displays the Nano-ITC control software interface. The top section is titled "Experiment Setup" and contains three radio buttons: "Incremental Titration" (selected), "Continuous Titration", and "Electrical Pulses". Below these are five buttons: "Insert..." with a green plus icon, "Delete" with a grey X icon, "Setup..." with a green and red icon, "Save As..." with a floppy disk icon, and "Load" with a folder icon. To the right of these buttons is a large table with three columns: "Injection", "Injection Interval (s)", and "Volume (μL)". The bottom section is titled "Equilibration" and contains two radio buttons: "Start Delay (s):" with a text box containing "0" and "(No data collected.)", and "Auto Equilibrate" (selected). Below "Auto Equilibrate" are three radio buttons for "Expected Heats": "Small", "Medium" (selected), and "Large". To the right of these is a checkbox for "Timeout (s):" which is checked, with a text box containing "1800".

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


The screenshot shows the 'Setup' window of the Nano-ITC software. The 'Run Mode Indicators' tab is selected. The 'Experiment Details' section (highlighted with a red box) contains the following fields:

- Syringe Concentration (mM): 1.3
- Cell Concentration (mM): 0.07
- Comments: RNAse into 2'CMP
- Sample Information**

The 'Experiment Equilibration Controls' section (highlighted with a red box) contains the following fields:

- Start Delay (s): 0 (No data collected)
- Auto Equilibrate: ☒ (Expected Heat: ☒ Small, ☐ Medium, ☐ Large)
- Timeout (s): 10
- Initial Baseline (s): 250 (Collected points before first injection)
- Final Baseline (s): 250 (Collected points after last injection)

The 'Experiment Event Timing' table (highlighted with a red box) shows the following data:

	Injection	Injection Interval (s)	Volume (μL)
	1	250	4
	2	250	4
	3	250	4
	4	250	4
	5	250	4
	6	250	4
	7	250	4
	8	250	4
	9	250	4
	10	250	4
	11	250	4
	12	250	4

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The screenshot shows the 'Experiment Equilibration Controls' dialog box in the Nano-ITC software. The dialog is titled 'Equilibration' and contains the following settings:

- Start Delay (s):** 0 (No data collected.)
- Auto Equilibrate:** Selected (radio button).
- Expected Heats:** Small (radio button), Medium (radio button), Large (radio button).
- Timeout (s):** 10
- Initial Baseline (s):** 250 (Collected points before first injection.)
- Final Baseline (s):** 250 (Collected points after last injection.)

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