


Type 100 Ti **BIOE**

8 x 6.0 mL

Type 100 Ti

**Fixed-Angle Rotor, Titanium**For use in instruments classified: **R S**

Major applications: Rapid separation of plasmid DNA and rapid differential centrifugation of small particles.

Max. RPM	Max. g	k Factor	Number of Tubes Volume/Size	Rotor Capacity
100 000	802 000	15	8 x 6 mL ½ x 2 ½ in 13 x 64 mm	48 mL

No. 363013. Type 100 Ti Rotor Assembly with 3 each 876089 O-rings, 3 each 839347 O-rings, and 2 each 363304 Overspeed Disks.

No. 363963. Quick-Seal Tube Kit with 4 boxes 344619 Polyallomer Tubes, 8 Spacers 365470, 1 each Tube Removal Tool 361668, Spinkote Lubricant 306812, and Vacuum Grease 335148. (Cordless Tube Topper and Tube Rack required.)

Tubes and Bottles

Tube Style/Material	Vol. (mL) g-Max	Part No.	Quantity	Size (mm)	Required Caps/ Spacers (pkg of 8)	Required Adapters (qty. 1)	g-Force	k Factor	Max. Speed
Quick-Seal Polyallomer	2.0	345829	50	13 x 25	360270	—	802 000	7.0	100 000
	2.0	363962 [§]	50	13 x 25	360270	—	802 000	7.0	100 000
	3.5	349621	50	13 x 32	360270	—	802 000	9.4	100 000
	3.5	363961 [§]	50	13 x 32	360270	—	802 000	9.4	100 000
	5.1	362248	50	13 x 51	360270	—	802 000	13.5	100 000
	5.1	363960 [§]	50	13 x 51	360270	—	802 000	13.5	100 000
	6.0	344619	50	13 x 64	365470	—	802 000	15.0	100 000

Note: To help you locate g-Max tubes in the rotor listings, these tubes are highlighted with bold, italic type.

[§] To simplify ordering, this g-Max Kit includes enough spacers to fill rotor, 50 tubes and required tools.

Rotor Replacement Parts

363037 Rotor Lid
 363039 Rotor Handle
 363304 Overspeed Disk, 100 000 rpm
 839347 O-ring for Rotor Handle
 876089 O-ring for Rotor Lid

Spacers (Package of 8)

365470 fits 13 x 64 mm
 360270 fits 13 x 51 mm
 fits 13 x 32 mm
 fits 13 x 25 mm

Adapters/Spacers/Caps

360270 365470

